







Vario Kugel SnapVKS-OC 1.7 mm

	43006550	01	(4 pounds of retained force)	Green
	43006590	02	(6 pounds of retained force)	Yellow
	43006560	03	(8 pounds of retained force)	Red
	430OC104 supersnap	04	(10 pounds of retained force)	Light Green
	430OC105 supersnap	05	(12 pounds of retained force)	Light Yellow
	430OC106 supersnap	06	(14 pounds of retained force)	Light Red
	430OC107 supersnap	07	(16 pounds of retained force)	Very Light Green




Vario Kugel SnapVKS-OC 2.2 mm

	43005440	01	(4 pounds of retained force)	Green
	43005450	02	(6 pounds of retained force)	Yellow
	43005460	03	(8 pounds of retained force)	Red
	430OC204 supersnap	04	(10 pounds of retained force)	Light Green
	430OC205 supersnap	05	(12 pounds of retained force)	Light Yellow
	430OC206 supersnap	06	(14 pounds of retained force)	Light Red




Bredent recommends for the following helpful hints to produce a successful long lasting prosthesis:

1. Attachment males must be cast in at least 270 Vickers Hardness to prevent premature wearing of the ball.
2. Never over-polish the male parts since this will hinder adequate retention.
3. For new cases alays begin with the lowest retention level first, which is Green Level # 1.
4. Remind your patients not to "bite" their prosthesis in place, but insert with finger pressure.
5. Patients need to clean their prosthesis at least twice daily; studies have shown that premature wear can occur to the prosthesis if not cleaned regularly.
6. Patients must see their dental health practitioner at least every six months to check the accuracy of the fit of their prosthesis, failing to do this will compromise the longevity of the case.
7. Should higher retention levels be needed, go up one level at a time to ensure patient comfort and ease of use.

Vario Kugel SnapVKS-SG 1.7 mm

	43006680	01	(4 pounds of retained force)	Green
	43006660	02	(6 pounds of retained force)	Yellow
	43006640	03	(8 pounds of retained force)	Red
	430SG104 supersnap	04	(10 pounds of retained force)	Light Green
	430SG105 supersnap	05	(12 pounds of retained force)	Light Yellow
	430SG106 supersnap	06	(14 pounds of retained force)	Light Red
	430SG107 supersnap	07	(16 pounds of retained force)	Very Light Green
	430SG109 supersnap	09	(20 pounds of retained force)	Very Light Red
	430SG110 supersnap	10	(22 pounds of retained force)	Lightest Green

Vario Kugel SnapVKS-SG 2.2 mm

	43005410	01	(4 pounds of retained force)	Green
	43005420	02	(6 pounds of retained force)	Yellow
	43005430	03	(8 pounds of retained force)	Red
	430SG204 supersnap	04	(10 pounds of retained force)	Light Green
	430SG205 supersnap	05	(12 pounds of retained force)	Light Yellow
	430SG206 supersnap	06	(14 pounds of retained force)	Light Red
	430SG208 supersnap	08	(18 pounds of retained force)	Very Light Yellow
	430SG209 supersnap	09	(20 pounds of retained force)	Very Light Red
	430SG210 supersnap	10	(22 pounds of retained force)	Lightest Green